

	ZOOM ® View™		Report Created 02 Jul 2014
	Device Settings Report		Last Office Interrogation N/R
	Date of Birth	N/R N/R N/R	Implant Date N/R
	Device	INVIVE W173/118030	

Programming

Last Programmed 02 Jul 2014

Ventricular Tachy

Ventricular Tachy EGM Storage On

Detection Rate 160 min⁻¹ (375 ms)

Atrial Tachy

Therapy

Ventricular Regulation

Vent Rate Regulation Off

BiV Trigger Off

Brady/CRT

Settings

Mode	VVIR
Lower Rate Limit	80 min ⁻¹
Maximum Sensor Rate	130 min ⁻¹
RV-Refractory (RVRP)	230 - 250 ms
LV-Refractory (LVRP)	250 ms
Ventricular Pacing Chamber	BiV
LV Offset	-40 ms
LV Protection Period	400 ms
Blanking	
A-Blank after V-Pace	125 ms
A-Blank after RV-Sense	45 ms
Magnet Response	Pace Async
Noise Response	VOO

Rate Enhancements

Rate Smoothing	
Up	Off %
Down	Off %

Output

■RV	Auto 2.0 V @ 0.4 ms
◆LV	3.0 V @ 0.4 ms

Sensitivity

●A	Fixed 0.75 mV
■RV	Fixed 2.5 mV
◆LV	Fixed 2.5 mV

Leads

●A	
Sense	Bipolar
■RV	
Pace	Bipolar
Sense	Bipolar
Safety Switch	On
◆LV	
Electrode Configuration	Dual
Pace	LVtip>>LVring
Sense	LVtip>>LVring
Safety Switch	On

Rate Adaptive Pacing

Accelerometer	On
Response Factor	8
Activity Threshold	Medium
Reaction Time	30 s
Recovery Time	2 min

Setup

Telemetry

Enable use of ZIP™ telemetry	On
------------------------------	----

Trending Setup

Recording Method	30 Second Average
Duration	25 hours
Data Storage	Continuous

Sleep Schedule

Sleep Start Time	23:00 hh:mm
Sleep Duration	07 hours

Leads Status Setup

	Intrinsic Amplitude	Pace Impedance	Pace Impedance Alert Range	Pace Threshold
●A	On	On	200 - 2000 Ω	
■RV	On	On	200 - 2000 Ω	On
◆LV	On	On	200 - 2000 Ω	

Other Daily Trends

Respiration-related Trends	On
----------------------------	----