

	ZOOM ® View™		Report Created 17 Mar 2017
	Device Settings Report		Last Office Interrogation 20 Oct 2016
	Date of Birth	N/R N/R N/R	Implant Date
	Device	VALITUDE X4 U128/705938	N/R

Programming

Last Programmed 17 Mar 2017

Ventricular Tachy

Ventricular Tachy EGM Storage On

Detection Rate 160 min⁻¹ (375 ms)

Atrial Tachy

Therapy

Ventricular Regulation

Vent Rate Regulation Off

BiV Trigger Off

Brady/CRT

Settings

Mode	VVIR	Output	
Lower Rate Limit	85 min ⁻¹	■RV	Trend 3.0 V @ 0.4 ms
Maximum Sensor Rate	130 min ⁻¹	◆LV	5.0 V @ 0.6 ms
RV-Refractory (RVRP)	230 - 250 ms	Sensitivity	
LV-Refractory (LVRP)	250 ms	●A	Fixed 0.75 mV
Ventricular Pacing Chamber	BiV	■RV	Fixed 2.5 mV
LV Offset	-30 ms	◆LV	Fixed 2.5 mV
LV Protection Period	400 ms	Leads	
Blanking		●A	
A-Blank after V-Pace	125 ms	Sense	Off
A-Blank after RV-Sense	45 ms	■RV	
Magnet Response	Pace Async	Pace	Bipolar
Noise Response	VOO	Sense	Bipolar
Rate Enhancements		Safety Switch	On
Rate Smoothing		◆LV	
Up	Off %	Electrode Configuration	Quadripolar
Down	Off %	Pace	LVTip1>>Can
		Sense	LVTip1>>LVRing2
		Rate Adaptive Pacing	
		Accelerometer	On
		Response Factor	8
		Activity Threshold	Medium
		Reaction Time	30 s
		Recovery Time	2 min

Setup

Telemetry

Enable use of ZIP™ telemetry On

Trending Setup

Recording Method 30 Second Average
Duration 25 hours
Data Storage Continuous

Sleep Schedule

Sleep Start Time 23:00 hh:mm
Sleep Duration 07 hours

Leads Status Setup

		Daily Intrinsic Amplitude	Daily Impedance	Impedance Limits Low High	Pace Threshold
●A	Pace/Sense	On	On	200 - 2000 Ω	Off
■RV	Pace/Sense	On	On	200 - 2000 Ω	On
◆LV	Pace/Sense	On	On	200 - 2000 Ω	

Other Daily Trends

Respiration-related Trends On