

	ZOOM ® View™		Report Created 07 Jul 2021
	Device Settings Report		Last Office Interrogation N/R
	Date of Birth	N/R N/R N/R	Implant Date N/R
	Device	ACCOLADE MRI L310/836133	

Programming

Last Programmed 07 Jul 2021

Tachy

Tachy EGM Storage On
Detection Rate 160 min⁻¹ (375 ms)

Atrial Tachy

Therapy

Ventricular Regulation

Vent Rate Regulation Off

Brady

Settings

Mode VVI
Lower Rate Limit 60 min⁻¹
Refractory 250 ms
Magnet Response Pace Async
Noise Response VOO

Rate Enhancements

Rate Smoothing
Down Off %
Rate Hysteresis
Hysteresis Offset Off min⁻¹

Output

Trend 3.0 V @ 0.4 ms

Sensitivity

Fixed 2.0 mV

Leads

Pace Bipolar
Sense Bipolar
Safety Switch On

Rate Adaptive Pacing

Minute Ventilation Passive
Accelerometer Passive

Minute Ventilation Sensor Settings

Minute Ventilation Sensor Passive
Excitation Current 320 µA
Vector Selection V Only
Signal Artifact Monitor On

Minute Ventilation Sensor Status

Status Calibrated
Current Vector V

Setup

Telemetry

Enable use of ZIP™ telemetry On

Trending Setup

Recording Method 30 Second Average
Duration 25 hours
Data Storage Continuous

Sleep Schedule

Sleep Start Time 23:00 hh:mm
Sleep Duration 07 hours

Leads Status Setup

	Daily Intrinsic Amplitude	Daily Impedance	Impedance Limits Low High	Pace Threshold
Pace/Sense	On	On	200 - 2000 Ω	On