

	ZOOM ® View™	Report Created 07 Jan 2019
	Device Settings Report Date of Birth N/R N/R N/R Device VITALIO MRI J275/101893	Last Office Interrogation 08 Jan 2018 Implant Date 2 Jun 2014

Programming

Last Programmed 07 Jan 2019

Tachy

Tachy EGM Storage On
 Detection Rate 160 min⁻¹ (375 ms)

Atrial Tachy

Therapy

Ventricular Regulation

Vent Rate Regulation Off

Brady

Settings

Mode	VVIR	Output	
Lower Rate Limit	55 min ⁻¹		Auto 1.3 V @ 0.4 ms
Maximum Sensor Rate	115 min ⁻¹	Sensitivity	
Refractory	250 - 250 ms		Fixed 2.0 mV
Magnet Response	Pace Async	Leads	
Noise Response	VOO	Pace	Unipolar
Rate Enhancements		Sense	Bipolar
Rate Smoothing		Rate Adaptive Pacing	
Up	Off %	Minute Ventilation	On
Down	Off %	Response Factor	8
		Fitness Level	Active
		Ventilatory Threshold	120 min ⁻¹
		Ventilatory Thresh. Response	70 %
		Accelerometer	Passive

Setup

Telemetry

Enable use of ZIP™ telemetry On

Trending Setup

Recording Method 30 Second Average
Duration 25 hours
Data Storage Continuous

Sleep Schedule

Sleep Start Time 23:00 hh:mm
Sleep Duration 07 hours

Leads Status Setup

	Daily Intrinsic Amplitude	Daily Impedance	Impedance Limits Low High	Pace Threshold
Pace/Sense	On	On	200 - 2000 Ω	On