Scientific Scientific

ZOOM ® View™		Report Created 21 Mar 2017
Combined Follo	ow-up Report	Last Office Interrogation
		13-JUL-2013
		Implant Date
Device	ALTRUA 50 S508/576134	N.R. N.R. N.R.

Events Since Last Reset

09-JUL-2013 12:36 Tach at 190 min⁻¹
09-JUL-2013 12:29 Tach at 280 min⁻¹
See last page for full list of Events Since Last Reset.

Battery Good

Intrinsic Amplitude

Pace Impedance

Pace Threshold

Longevity Remaining >5.0 years

@ Current Pacing Percentage(s)

Magnet Rate 100 min⁻¹

Elective Replacement (ERT)
End of Life (EOL)

End of Life (EOL) (BOL)

Previous Session Present Session

13-JUL-2013

3.3 mV N.R. mV 560Ω 490 Ω 1.3 V @ 0.6 ms

Beginning of Life

Settings

Leads Data

Tachy Settings

Tachy Detection On
Rate 170 min⁻¹
Duration 8 cycles

Brady Settings Pacing Output 0.1 V @ 0.05 ms Mode 2.0 mV SSI Sensitivity Lead Configuration Lower Rate Limit 30 min⁻¹ **Bipolar** -- min⁻¹ Pace Maximum Tracking Rate Refractory 250 ms Sense

Implant

N.R. N.R. N.R.


N.R. mV

N.R. Ω


N.R. V

Daily Measurement

Date	Amplitude	Impedance	Threshold
	(mV)	(Ω)	(V)
12-JUL-2013	N.R.	490	Off
11-JUL-2013	N.R.	480	Off
10-JUL-2013	N.R.	520	Off
09-JUL-2013	N.R.	550	Off
08-JUL-2013	3.5	600	Off
07-JUL-2013	8.0	640	Off
06-JUL-2013	N.R.	610	Off
04-JUL-2013	3.5	590	Off
27-JUN-2013	8.3	600	Off
20-JUN-2013	7.6	610	Off
13-JUN-2013	8.1	620	Off
06-JUN-2013	8.0	610	Off
30-MAY-2013	5.6	530	Off



Tachy Counters			Since Last Reset
Episode Counters			
Tachy Detections			5
Brady Counters		Percentages	Since Last Reset
Counters			
Paced		87 %	3.4M
Sensed		13 %	499.0K
Intrinsic Promotion			
Rate Hysteresis			
Searches			0
Successful Searches			0
Counters			
Interval Variation			
Variation 0 <= 10%			3.9M
Variation 10 <= 20%			17.5K
Variation 20 <= 30%			3.9K
Variation > 30%			7.2K
Histograms	Since Last Reset		



All Events Since Last Reset

09-JUL-2013 12:36	Tach at 190 min ⁻¹
09-JUL-2013 12:29	Tach at 280 min ⁻¹
09-JUL-2013 12:25	Tach at 240 min ⁻¹
09-JUL-2013 12:25	Tach at 260 min ⁻¹
09-JUL-2013 12:23	Tach at 220 min ⁻¹

2892 Software Version: 1.00 S508 Firmware Version: 0a © 2010

Boston Scientific Corporation or its affiliates. All rights reserved.

Page 3 of 3

Clinician Signature: