

ZOOM ® View™ Report Created 22 Jul 2015
---

### **Combined Follow-up Report**

Date of Birth N/R N/R N/R

Device AUTOGEN ICD D174/191431

Tachy Mode Monitor + Therapy

Last Office Interrogation

N/R

Implant Date N/R

#### My Alerts

There are no alerts to display.

#### **Events Since Last Reset (22 Jul 2015)**

There are no events to display.

Battery OK

Approximate time to explant: > 8 years
Charge Time N/R s

Last Capacitor Re-form N/R

One Year Remaining

**Explant Leads Data Implant Previous Session Most Recent** N/R N/R N/R Ventricular Intrinsic Amplitude N/R mV N/R mV 15.2 mV @ 71 min<sup>-1</sup> Pace Impedance  $N/R \Omega$  $N/R \Omega$ 526 Ω Pace Threshold N/R V @ N/R ms N/R Auto 0.4 V @ 0.4 ms **Shock Vector** Shock Impedance  $N/R \Omega$  $N/R \Omega$ 74 Ω

Settings

**Ventricular Tachy Settings** 

VF 230 min<sup>-1</sup> ATP 41J, 41J, 41Jx6 VT 160 min<sup>-1</sup> Burst Ramp 41J, 41J, 41Jx4

VT-1 140 min<sup>-1</sup> Monitor Only

**Brady Settings** 

Mode VVI Pacing Output Lower Rate Limit 40 min<sup>-1</sup> Ventricular

V-Refractory (VRP) 250 ms Sen

Sensitivity
Ventricular
AGC 0.6 mV

Leads Configuration (Pace/Sense)

Ventricular Bipolar

3.5 V @ 0.4 ms

Rate Adaptive Pacing

Minute Ventilation Passive
Accelerometer Passive

2868 Software Version: 3.05

D174 Firmware Version: E\_v1.02.00(1.01)

© 2014

Boston Scientific Corporation
or its affiliates. All rights reserved.

Page 1 of 4

Clinician Signature:

Rhythm ID Update					
Last Successful Con	npletion		22 Jul 2015 11:38		
Trends					
Events		Apr 15	May 15	Jun 15	Jul 15
	V				
	Other				
Activity Level		Apr 15	May 15	Jun 15	Jul 15
•					
	%				
		A15			lul ar
Heart Rate Max		Apr 15	May 15	Jun 15	Jul 15
Mean					
Min	min <sup>-1</sup>				
Respiratory Rate		Apr 15	May 15	Jun 15	Jul 15
Max					
Median	rpm				
Min					
			! 		
AP Scan		Apr 15	May 15	Jun 15	Jul 15
	events/h				
R-Wave (V)		Apr 15	May 15	Jun 15	Jul 15
Amplitude <sup>*</sup>					
	mV				
V Pace Impedance		Apr 15	May 15	Jun 15	Jul 15
	•				
	Ω				
V Pace Threshold		Apr 15	May 15	Jun 15	Jul 15
	V				
Chaald lynnadan -		Apr 15	May 15	lue 15	iul 15
Shock Impedance		∨hı 13	May 15	Jun 15	Jul 15
	Ω				
	77				

2868 Software Version: 3.05 D174 Firmware Version: E\_v1.02.00(1.01) © 2014


Boston Scientific Corporation
or its affiliates. All rights reserved.

Page 2 of 4

Clinician Signature:

Ventricular Tachy Counters	Since Last Reset 22 Jul 2015	Device Totals	
Ventricular Episode Counters			
Total Episodes	0	0	
Treated			
VF Therapy	0	0	
VT Therapy	0	0	
VT-1 Therapy	0	0	
Commanded Therapy	0	0	
Non-Treated			
No Therapy Programmed	0	0	
Nonsustained Episodes	0	0	
Other Untreated Episodes	0	0	
Ventricular Therapy Counters			
ATP Delivered	0	0	
ATP % Successful	0	0	
Shocks Delivered	0	0	
First Shock % Successful	0	0	
Shocks Diverted	0	0	
Brady Counters	Reset Before Last	Since Last Reset	
	N/R	22 Jul 2015	
Counters			
% V Paced	0	0	
Intrinsic Promotion			
Rate Hysteresis			
% Successful	0	0	

## **Histograms Reset Before Last (N/R)**



# **Histograms Since Last Reset (22 Jul 2015)**

