



ZOOM ® View™

Device Settings Report

HOVADÍK 510404, Bronislav

Date of Birth 4 Apr 1951
 Device CHARISMA EL ICD D332/243151
 Tachy Mode Monitor + Therapy

Report Created 12 Nov 2018

Last Office Interrogation
08 Nov 2018
 Implant Date
8 Nov 2018

Programming

Last Programmed 12 Nov 2018
 Ventricular Tachy Mode 08 Nov 2018 Changed to Monitor + Therapy
 08 Nov 2018 Changed to Off

Ventricular Tachy

VF 220 min⁻¹ (273 ms)

Detection/Redetection

Initial Duration 2.5 s
 Redetection Dur 1.0 s
 Post-shock Dur 1.0 s

Therapy

QUICK CONVERT™ ATP 250 min⁻¹
 Shock 1 36 J
 Shock 2 41 J
 Additional 41 J Shocks 6

VT 170 min⁻¹ (353 ms)

Detection/Redetection

Initial Duration 2.5 s
 Redetection Dur 3.0 s
 Post-shock Dur 1.0 s

ATP1

Burst 2
 Number of Bursts
 Pulses per Burst
 Initial 10
 Increment 0
 Coupling Interval 81 %
 Decrement 0 ms
 Burst Cycle Length 81 %
 Ramp Decrement 0 ms
 Scan Decrement 0 ms
 Minimum Interval 280 ms

Enhancements

VT Detection Rhythm ID On
 Initial Detection On
 Sustained Rate Duration Off mm:ss

Post-Shock Detection

Rhythm ID Setup

Passive Method On
 Active Method On
 Temporary LRL 50 min⁻¹

Common Parameters

Stability (For Post-Shock only) 30 ms
 RhythmMatch™ Threshold 94 %

ATP2

Ramp/Scan
 Number of Bursts 2
 Pulses per Burst
 Initial 10
 Increment 0
 Coupling Interval 81 %
 Decrement 10 ms
 Burst Cycle Length 81 %
 Ramp Decrement 10 ms
 Scan Decrement 8 ms
 Minimum Interval 250 ms
 ATP Time-out Off mm:ss

Shocks

Shock 1 26 J
 Shock 2 31 J
 Shock 3 -6 41 J

Ventricular Tachy (Continued)

VT-1 150 min⁻¹ (400 ms)

Detection/Redetection

Initial Duration 60.0 s
Redetection Dur 3.0 s
Post-shock Dur 1.0 s

Enhancements

Rhythm ID

ATP1

Number of Bursts

Off

Off

ATP2

Number of Bursts

Off

Off

Shocks

Shock 1

Off J

Shock 2

Off J

Shock 3 -5

Off J

Ventricular Tachy Therapy Setup

ATP

Ventricular ATP Amplitude 5.0 V
Ventricular ATP Pulse Width 1.0 ms

Magnet and Beeper

Magnet Response Inhibit Therapy
Beep During Capacitor Charge Off

Shock (All Shocks)

Waveform Biphasic

Committed Shock Off

Lead Polarity Initial

Shock Lead Vector RV Coil to RA Coil
and Can

Atrial Tachy

Therapy

VTR Fallback LRL 70 min⁻¹

Ventricular Regulation

Vent Rate Regulation Off

Brady

Normal Settings

Mode	VVI	Output	■V	2.6 V @ 0.4 ms
Lower Rate Limit	50 min ⁻¹	Sensitivity	■V	AGC 0.6 mV
V-Refractory (VRP)	250 ms	Leads	■V	
Noise Response	VOO			
Rate Enhancements				
Rate Smoothing				
Up	Off %	Pace		Bipolar
Down	Off %	Sense		Bipolar
Rate Hysteresis		Rate Adaptive Pacing		
Hysteresis Offset	-10 min ⁻¹	Accelerometer		Passive
Search Hysteresis	Off cycles			

Brady (Post-Therapy)

Brady Settings

Lower Rate Limit	60 min ⁻¹	Post Therapy	
Output		Post Therapy Period	00:30 mm:ss
■V	5.0 V @ 1.0 ms		

Setup

Beeper

Beep when Explant is Indicated	On
--------------------------------	----

Telemetry

Enable use of ZIP™ telemetry	On
------------------------------	----

Trending Setup

Recording Method	30 Second Average
Duration	25 hours
Data Storage	Continuous

Sensor Activation

Heart Failure Sensor Suite	On
Sleep Incline	On

Sleep Schedule

Sleep Start Time	23:00 hh:mm
Sleep Duration	07 hours

Leads Status Setup

	Daily Intrinsic Amplitude	Daily Impedance	Impedance Limits Low High	Beep When Out-of-Range	Pace Threshold
■V Pace/Sense	On	On	200 - 2000 Ω	Off	Off
Shock		On	20 - 125 Ω	Off	

Other Daily Trends

Respiration-related Trends	On
----------------------------	----