

Combined Follow-up Report

Date of Birth N/R N/R N/R

Device ACCOLADE MRI L311/269927

Last Office Interrogation

N/R Implant Date N/R

My Alerts

There are no alerts to display.

Events Since Last Reset (23 Mar 2017)

26 Mar 2017 19:10 RA Auto 26 Mar 2017 18:16 RV Auto

Battery OK

Approximate time to explant: > 8 years
Magnet Rate 100 min⁻¹

One Year Remaining

Explant

		=×p.u	-
Leads Data	Implant N/R N/R N/R	Previous Session	Most Recent
Atrial			
Intrinsic Amplitude	N/R mV	N/R mV	2.2 mV
Pace Impedance	N/R Ω	N/R Ω	615 Ω
Pace Threshold	N/R V @ N/R ms	N/R	Auto 1.7 V @ 0.4 ms
Ventricular			
Intrinsic Amplitude	N/R mV	N/R mV	14.8 mV @ 40 min ⁻¹
Pace Impedance	N/R Ω	N/R Ω	687 Ω
Pace Threshold	N/R V @ N/R ms	N/R	Auto 0.7 V @ 0.4 ms
Cattings			

Settings			
1/	 T -	- 1	0 - 11

Detection Rate

Ventricular Tachy SettingsAtrial Tachy SettingsVentricular Tachy EGM StorageOnATR Mode Switch 170 min⁻¹ DDI

160 min⁻¹

Brady Settings

Mode	DDDR	Pacing Output	
RYTHMIQ™	Off	Atrial	Trend 3.0 V @ 0.4 ms
Lower Rate Limit	55 min ⁻¹	Ventricular	Trend 3.0 V @ 0.4 ms
Maximum Tracking Rate	130 min ⁻¹	Sensitivity	
Maximum Sensor Rate	130 min ⁻¹	Atrial	Fixed 0.75 mV
Paced AV Delay	80 - 180 ms	Ventricular	Fixed 2.5 mV
Sensed AV Delay	65 - 150 ms	Leads Configuration (Pace/S	Sense)
A-Refractory (PVARP)	240 - 280 ms	Atrial	Bipolar
V-Refractory (VRP)	230 - 250 ms	Ventricular	Bipolar
		Rate Adaptive Pacing	

Minute Ventilation On
Accelerometer Passive

2869 Software Version: 2.02 L311 Firmware Version: F_v1.00.00 © 2014

Boston Scientific Corporation
or its affiliates. All rights reserved.

Page 1 of 7

Clinician Signature:

Trends				
Today's Activity: F	Program and			
R = Counter Reset P = Program I = Interrogate		Jan 17	, Feb 17 ,	Mar 17
Events	A V Other		Feb 17	Mar 17
Activity Level	5 3.8 % 2.5 1.2		Feb 17	Mar 17
Heart Rate - Max - Mean - Min	min ⁻¹	Jan 17	Feb 17	Mar 17
SDANN	ms	Jan 17	Feb 17	Mar 17
HRV Footprint	%	Jan 17	Feb 17	Mar 17
Autonomic Balance Monitor		Jan 17	Feb 17	Mar 17
AT/AF Burden	4 3 h 2 1		Feb 17	Mar 17
V Rate during AT/AF Max Mean	min ⁻¹	Jan 17	Feb 17	Mar 17

2869 Software Version: 2.02 L311 Firmware Version: F_v1.00.00 © 2014

Boston Scientific Corporation
or its affiliates. All rights reserved.

Page 2 of 7

Clinician Signature:	

Trends (Continued)	1			
Pacing Percent - Atrial - V	100 95 % 90 85 80		Feb 17	15
Respiratory Rate - Max - Median - Min	32.5		Feb 17	
AP Scan	30 22.5 events/h 15 7.5	Jan 17	Feb 17	\\\

P-Wave Amplitude			Jan 17	Feb 17	Mar 17
	тV	3 2.8 2.5 2.2			
A Pace Impedance		21	Jan 17	Feb 17	Mar 17
	Ω	700 675 650 625 600	J 2.		Λ
A Pace Threshold		2.5	Jan 17	Feb 17	Mar 17
	V	2.5 1.9 1.2 0.6 0			
R-Wave (V)			Jan 17	Feb 17	Mar 17
Amplitude		22 21.5			•
	тV	21			
		20.5			
/ Pace Impedance		000	Jan 17	Feb 17	Mar 17
		800 750			
	Ω	700			```
		650 600			
/ Pace Threshold		2.5	Jan 17	Feb 17	Mar 17
		1.9			
	V	1.2			
		0.6			

Ventricular Tachy Counters	Since Last Reset 23 Mar 2017	Device Totals
Ventricular Episode Counters		
Total Episodes	0	0
VT Episodes (V>A)	0	0
SVT Episodes (V<=A)	0	0
Nonsustained Episodes	0	0
Brady Counters	Reset Before Last N/R day(s) N/R to 23 Mar 2017	Since Last Reset 4 day(s) 23 Mar 2017 to Today
Counters		
% A Paced	0	96
% V Paced	0	100
Intrinsic Promotion		
AV Search +		
% Successful	0	0
Rate Hysteresis		
% Successful	0	0
Atrial Arrhythmia		
% AT/AF	N/R	0
Total Time in AT/AF (min)	N/R	0.0
Episodes by Duration		
< 1 minute	0	0
1 min - < 1 hr	0	0
1 hr - < 24 hr	0	0
24 hr - < 48 hr	0	0
> 48 hr	0	0
Total PACs	0	710
Ventricular Counters		
Total PVCs	0	774
Three or More PVCs	0	3
AT/AF Overview: Since Last Reset (23 Mar 201	7 to Today)	

AT/AF Overview: Since Last Reset (23 Mar 2017 to Today)

AT/AF: 0% AT/AF Events: 0 Total Time in AT/AF (min): 0.0 ATR Episodes initiated within 2 hours of lead attach are not stored in Arrhythmia Logbook.


AT/AF Overview: Reset Before Last (N/R to 23 Mar 2017)


ATR Episodes initiated within 2 hours of lead attach are not stored in Arrhythmia Logbook.


2869 Software Version: 2.02 L311 Firmware Version: F_v1.00.00

Clinician Signature:


Histograms Reset Before Last (N/R)







Histograms Since Last Reset (23 Mar 2017)




2869 Software Version: 2.02 L311 Firmware Version: F_v1.00.00 © 2014


Boston Scientific Corporation
or its affiliates. All rights reserved.

Page 6 of 7

Clinician Signature:

Histograms Since Last Reset (23 Mar 2017) (Continued)





Heart Rate Variability

Last Measured:

HRV data unavailable because device was reprogrammed in the last 24 hours.

Reference:

HRV data unavailable because device was reprogrammed in the last 24 hours.