



ZOOM® View™

Device Settings Report

KRÍZ, Karel

Date of Birth 30 Aug 1954
 Device AUTOGEN X4 CRT-D G179/136339
 Tachy Mode Monitor + Therapy

Report Created 27 Jun 2016

Last Office Interrogation
03 Jun 2016
 Implant Date
12 Oct 2015

Programming

Last Programmed 27 Jun 2016
 Ventricular Tachy Mode 12 Oct 2015 Changed to Monitor + Therapy
 12 Oct 2015 Changed to Off

Ventricular Tachy

VF 230 min⁻¹ (261 ms)

Detection/Redetection

Initial Duration 2.0 s
 Redetection Dur 1.0 s
 Post-shock Dur 1.0 s

Therapy

QUICK CONVERT™ ATP 250 min⁻¹
 Shock 1 36 J
 Shock 2 41 J
 Additional 41 J Shocks 6

VT 185 min⁻¹ (324 ms)

Detection/Redetection

Initial Duration 4.0 s
 Redetection Dur 1.0 s
 Post-shock Dur 1.0 s

ATP1

Burst 2
 Number of Bursts 2
 Pulses per Burst Initial 10
 Increment 0
 Coupling Interval 81 %
 Decrement 0 ms
 Burst Cycle Length 81 %
 Ramp Decrement 0 ms
 Scan Decrement 0 ms
 Minimum Interval 220 ms

Enhancements

Rhythm ID
 VT Detection On
 Initial Detection On
 Sustained Rate Duration 05:00 mm:ss
 Post-Shock Detection On
 Sustained Rate Duration 02:00 mm:ss

Rhythm ID Setup

Ambulatory update On
 Temporary LRL 85 min⁻¹

Common Parameters

Atrial Tachy Discrimination On
 AFib Rate Threshold 170 min⁻¹
 Stability 20 ms
 RhythmMatch™ Threshold 94 %

ATP2

Ramp 2
 Number of Bursts 2
 Pulses per Burst Initial 10
 Increment 0
 Coupling Interval 78 %
 Decrement 0 ms
 Burst Cycle Length 78 %
 Ramp Decrement 10 ms
 Scan Decrement 0 ms
 Minimum Interval 220 ms
 ATP Time-out 01:00 mm:ss

Shocks

Shock 1 31 J
 Shock 2 36 J
 Shock 3 -6 41 J

Ventricular Tachy Therapy Setup

ATP

RV ATP Amplitude 5.0 V
 RV ATP Pulse Width 1.0 ms
 LV ATP Amplitude 5.0 V
 LV ATP Pulse Width 1.0 ms

Shock (All Shocks)

Waveform Biphasic
 Committed Shock Off
 Lead Polarity Initial
 Shock Lead Vector RV Coil to RA Coil and Can

Magnet and Beeper

Magnet Response Inhibit Therapy

Ventricular Tachy (Continued)

Ventricular Tachy Therapy Setup

Magnet and Beeper

Beep During Capacitor Charge Off

Atrial Tachy

Therapy

ATR Mode Switch Details

ATR Mode Switch	On
Trigger Rate	170 min ⁻¹
Duration	8 cycles
Entry Count	8 cycles
Exit Count	8 cycles
Fallback	
Mode	DDI
Time	00:30 mm:ss
ATR/VTR Fallback LRL	70 min ⁻¹
Ventricular Rate Regulation	Min
BiV Trigger	On
Maximum Pacing Rate	130 min ⁻¹

Atrial Tachy Response

Atrial Flutter Response	Off
PMT Termination	On

Ventricular Regulation

BiV Trigger	On
-------------	----

Brady/CRT

Normal Settings

Mode	DDDR	Output	
Lower Rate Limit	85 min ⁻¹	●A	2.5 V @ 0.4 ms
Maximum Tracking Rate	130 min ⁻¹	■RV	3.0 V @ 0.4 ms
Maximum Sensor Rate	130 min ⁻¹	◆LV	3.5 V @ 0.4 ms
Paced AV Delay	220 - 250 ms	Sensitivity	
Sensed AV Delay	165 - 190 ms	●A	AGC 0.25 mV
A-Refractory (PVARP)	240 - 280 ms	■RV	AGC 0.6 mV
RV-Refractory (RVRP)	230 - 250 ms	◆LV	AGC 1.0 mV
LV-Refractory (LVRP)	250 ms	Leads	
Ventricular Pacing Chamber	BiV	●A	
LV Offset	-60 ms	Pace	Bipolar
PVARP after PVC	400 ms	Sense	Bipolar
LV Protection Period	400 ms	■RV	
Blanking		Pace	Bipolar
A-Blank after V-Pace	Smart ms	Sense	Bipolar
A-Blank after RV-Sense	Smart ms	◆LV	
RV-Blank after A-Pace	65 ms	Electrode Configuration	Quadripolar
LV-Blank after A-Pace	Smart ms	Pace	LVRing3>>Can
Noise Response	DOO	Sense	LVTip1>>LVRing2
Rate Enhancements		Rate Adaptive Pacing	
Rate Smoothing		Minute Ventilation	On
Up	Off %	Response Factor	8
Down	9 %	Fitness Level	Active
Tracking Preference	On	Ventilatory Threshold	125 min ⁻¹
		Ventilatory Thresh. Response	70 %
		Accelerometer	Passive

Brady/CRT (Post-Therapy)

Brady/CRT Settings

Lower Rate Limit	45 min ⁻¹	Post Therapy	
Output		Post Therapy Period	00:30 mm:ss
●A	5.0 V @ 1.0 ms		
■RV	5.0 V @ 1.0 ms		
◆LV	5.0 V @ 1.0 ms		

Setup

Beeper

Beep when Explant is Indicated On

Telemetry

Enable use of ZIP™ telemetry On

Trending Setup

Recording Method 30 Second Average
Duration 25 hours
Data Storage Continuous

Sleep Schedule

Sleep Start Time 23:00 hh:mm
Sleep Duration 07 hours

Leads Status Setup

	Daily Intrinsic Amplitude	Daily Impedance	Impedance Limits Low High	Beep When Out-of-Range	Pace Threshold
●A Pace/Sense	On	On	200 - 2000 Ω	Off	Off

Setup (Continued)

Leads Status Setup

		Daily Intrinsic Amplitude	Daily Impedance	Impedance Limits		Beep When Out-of-Range	Pace Threshold
				Low	High		
■RV	Pace/Sense	On	On	200	2000 Ω	Off	Off
◆LV	Pace/Sense	On	On	200	2000 Ω	Off	Off
	Shock		On	20	125 Ω	Off	