

ZOOM ® View™	Report Created 06 Jun 2016

**Combined Follow-up Report** 

Hlucík, Stepán Date of Birth

Device

7 Jan 1944 INCEPTA CRT-D P162/119874

Monitor + Therapy Tachy Mode

Last Office Interrogation

12 Oct 2015 Implant Date 12 Mar 2013

My Alerts

**a** 06 Jun 2016 13:55 Check Atrial and LV Leads

**Events Since Last Reset (14 Mar 2013)** 

VT at 200 min<sup>-1</sup>, ATPx4, 17J 30 Aug 2014 13:14 VT at 181 min<sup>-1</sup>, **♦** ATPx4, 17J 09 Aug 2014 16:46

See last page for full list of Events Since Last Reset.

OK **Battery** 

> Approximate time to explant: 5 years Charge Time 9.9 s

> Last Capacitor Re-form 26 Apr 2016 09:20



		Explain.	•
Leads Data	<b>Implant</b> 12 Mar 2013	Previous Session	Most Recent
Atrial			
Intrinsic Amplitude	N/R mV	N/R mV	N/R mV
Pace Impedance	N/R $\Omega$	>2500 Ω	>2500 Ω
Pace Threshold	N/R V @ N/R ms	4.0 V @ 0.4 ms	
Right Ventricular			
Intrinsic Amplitude	7.0 mV	8.1 mV	8.0 mV @ 71 min <sup>-1</sup>
Pace Impedance	$650~\Omega$	518 Ω	556 $\Omega$
Pace Threshold	0.5 V @ 0.4 ms	1.1 V @ 0.4 ms	1.3 V @ 0.4 ms
Left Ventricular			
Intrinsic Amplitude	12.0 mV	13.7 mV	19.0 mV
Pace Impedance	$650~\Omega$	512 Ω	$532~\Omega$
Pace Threshold	1.3 V @ 0.4 ms	1.3 V @ 0.4 ms	1.4 V @ 0.4 ms
Shock Vector			
Shock Impedance	45 Ω	55 Ω	57 Ω
Settings			

Ventricular Tachy Settings		
Vontrioular Tachy Sottings		

VF	230 min <sup>-1</sup> ATP		31J, 41J, 41Jx6
VT	205 min⁻¹ Scan	Ramp	31J, 41J, 41Jx4
VT-1	1 160 min <sup>-1</sup> Scan	Ramp	41J. 41J. 41Jx3

#### **Brady/CRT Settings**

Mode	VVIR	Pacing Output	
Lower Rate Limit	60 min <sup>-1</sup>	Right Ventricular	3.0 V @ 0.4 ms
Maximum Sensor Rate	130 min <sup>-1</sup>	Left Ventricular	3.0 V @ 0.4 ms
RV-Refractory (RVRP)	230 - 250 ms	Sensitivity	
LV-Refractory (LVRP)	250 ms	Atrial	AGC 0.25 mV
Ventricular Pacing Chamber	BiV	Right Ventricular	AGC 0.6 mV
LV Offset	0 ms	Left Ventricular	AGC 1.0 mV

Leads Configuration (Pace/Sense)

Bipolar/Off Atrial Right Ventricular **Bipolar** LVring>>RV / LVtip>>LVring Left Ventricular


Sensors


Accelerometer On

2868 Software Version: 3.07 P162 Firmware Version: B\_v1.02.00(4.01)

© 2014 Boston Scientific Corporation or its affiliates. All rights reserved. Page 1 of 6


Clinician Signature:







Ventricular Tachy Counters	Since Last Reset 14 Mar 2013	Device Totals	
Ventricular Episode Counters			
Total Episodes	7.1K	7.1K	
Treated			
VF Therapy	0	0	
VT Therapy	4	4	
VT-1 Therapy	17	17	
Commanded Therapy	0	0	
Non-Treated			
No Therapy Programmed	0	0	
Nonsustained Episodes	7.1K	7.1K	
Other Untreated Episodes	7	7	
Ventricular Therapy Counters			
ATP Delivered	29	29	
ATP % Successful	67	67	
Shocks Delivered	6	6	
First Shock % Successful	100	100	
Shocks Diverted	2	2	
Brady/CRT Counters	Reset Before Last 12 Mar 2013	Since Last Reset 14 Mar 2013	
Counters			
% A Paced	100	0	
% RV Paced	90	91	
% LV Paced	99	98	
Intrinsic Promotion			
Rate Hysteresis			
% Successful	0	0	
Atrial Burden			
Atrial Burden %	N/R	0	
(One Year of Data)			
Episodes by Duration			
< 1 minute	0	0	
1 min - < 1 hr	0	0	
1 hr - < 24 hr	0	0	
24 hr - < 48 hr	0	0	
> 48 hr	0	0	
Total PACs	0	0	
Ventricular Counters			
Total PVCs	14	69	
Three or More PVCs	1	6	


# **Histograms Reset Before Last (12 Mar 2013)**





## **Histograms Since Last Reset (14 Mar 2013)**





2868 Software Version: 3.07 P162 Firmware Version: B\_v1.02.00(4.01) © 2014

Boston Scientific Corporation
or its affiliates. All rights reserved.

Page 5 of 6

Clinician Signature:	

# **Heart Rate Variability**

### Last Measured:

Last Measured Data unavailable

### Reference:

Reference Data Unavailable

## All Events Since Last Reset (14 Mar 2013)

Events Since Last Reset	(14 Mar 2013)
03 Jun 2016 16:58	NonSustV at 159 min <sup>-1</sup> , Nonsustained
03 Jun 2016 15:21	NonSustV at 155 min <sup>-1</sup> , Nonsustained
03 Jun 2016 15:19	NonSustV at 153 min <sup>-1</sup> , Nonsustained
03 Jun 2016 15:15	NonSustV at 157 min <sup>-1</sup> , Nonsustained
03 Jun 2016 15:15	NonSustV at 156 min <sup>-1</sup> , Nonsustained
03 Jun 2016 15:14	NonSustV at 159 min <sup>-1</sup> , Nonsustained
03 Jun 2016 15:07	NonSustV at 161 min <sup>-1</sup> , Nonsustained
03 Jun 2016 15:02	NonSustV at 160 min <sup>-1</sup> , Nonsustained
03 Jun 2016 15:00	NonSustV at 155 min <sup>-1</sup> , Nonsustained
03 Jun 2016 14:59	NonSustV at 159 min <sup>-1</sup> , Nonsustained
03 Dec 2015 12:23	VT-1 at 181 min <sup>-1</sup> , No Therapy
03 Dec 2015 12:11	VT-1 at 172 min <sup>-1</sup> , No Therapy
03 Dec 2015 11:32	VT-1 at 146 min <sup>-1</sup> , No Therapy
19 Oct 2015 14:51	VT-1 at 175 min <sup>-1</sup> , No Therapy
23 Sep 2015 16:09	VT-1 at 163 min <sup>-1</sup> , ATPx1
23 Sep 2015 15:54	VT-1 at 173 min <sup>-1</sup> , ATPx1
23 Sep 2015 14:11	VT-1 at 157 min <sup>-1</sup> , ATPx1
23 Sep 2015 13:11	VT-1 at 162 min <sup>-1</sup> , <b>♦</b> ATPx4
23 Sep 2015 13:10	VT-1 at 179 min⁻¹, <b>◈</b> ATPx5
23 Sep 2015 13:07	VT-1 at 173 min <sup>-1</sup> , <b></b> ATPx2
21 Sep 2015 20:26	VT-1 at 178 min <sup>-1</sup> , ATPx1
21 Sep 2015 19:19	VT-1 at 186 min <sup>-1</sup> , ATPx1
17 Sep 2015 12:20	VT-1 at 181 min <sup>-1</sup> , <b>♦</b> ATPx8, 41J
17 Sep 2015 12:18	VT-1 at 181 min <sup>-1</sup> , <b></b> ATPx5
17 Sep 2015 10:44	VT-1 at 196 min <sup>-1</sup> , ATPx1
30 Aug 2015 17:02	VT-1 at 155 min <sup>-1</sup> , <b></b> ATPx3
30 Aug 2015 16:19	VT-1 at 165 min <sup>-1</sup> , No Therapy
29 Aug 2015 14:33	VT-1 at 170 min <sup>-1</sup> , No Therapy
29 Aug 2015 14:32	VT-1 at 186 min <sup>-1</sup> , ATPx1
29 Aug 2015 14:20	VT-1 at 176 min <sup>-1</sup> , ATPx1
18 May 2015 10:09	VT-1 at 180 min <sup>-1</sup> , No Therapy
25 Oct 2014 11:49	VT-1 at 183 min <sup>-1</sup> , <b>♦</b> ATPx8, 41J
25 Oct 2014 11:38	VT-1 at 155 min <sup>-1</sup> , ATPx1
10 Oct 2014 12:13	VT-1 at 168 min <sup>-1</sup> , <b>♦</b> ATPx8, 41J
30 Aug 2014 13:14	VT at 200 min⁻¹, <b>◈</b> ATPx4, 17J
09 Aug 2014 16:46	VT at 181 min⁻¹, <b></b> ATPx4, 17J
16 May 2013 17:54	VT at 178 min⁻¹, <b>ጭ</b> ATPx2
17 Apr 2013 15:14	VT at 179 min⁻¹, <b>ጭ</b> ATPx4, 17J

2868 Software Version: 3.07 P162 Firmware Version: B\_v1.02.00(4.01) © 2014
Boston Scientific Corporation
or its affiliates. All rights reserved.
Page 6 of 6

Clinician Signature: