



ZOOM® View™

Device Settings Report

Report Created 02 Jun 2014

Date of Birth N/R N/R N/R
Device VITALIO MRI J276/102818

Last Office Interrogation
N/R
Implant Date
N/R

Programming

Last Programmed 02 Jun 2014

Ventricular Tachy

Ventricular Tachy EGM Storage On
Detection Rate 160 min⁻¹ (375 ms)

Atrial Tachy

Therapy

ATR Mode Switch Details

ATR Mode Switch On
Trigger Rate 170 min⁻¹
Duration 8 cycles
Entry Count 8 cycles
Exit Count 8 cycles
Fallback
Mode DDI
Time 00:30 mm:ss
ATR Fallback LRL 70 min⁻¹
Ventricular Rate Regulation On
Maximum Pacing Rate 130 min⁻¹

Atrial Tachy Response

Atrial Flutter Response On
Trigger Rate 170 min⁻¹
PMT Termination On
Atrial Overdrive
APP/ProACt Off

Brady

Settings

| | | | |
|--------------------------|-----------------------|-----------------------------|----------------|
| Mode | DDD | Output | |
| RYTHMIQ™ | Off | ●A | 2.5 V @ 0.4 ms |
| Lower Rate Limit | 55 min ⁻¹ | ■V | 2.5 V @ 0.4 ms |
| Maximum Tracking Rate | 130 min ⁻¹ | Sensitivity | |
| Paced AV Delay | 80 - 180 ms | ●A | Fixed 0.75 mV |
| Sensed AV Delay | 65 - 150 ms | ■V | Fixed 2.5 mV |
| A-Refractory (PVARP) | 240 - 280 ms | Leads | |
| V-Refractory (VRP) | 230 - 250 ms | ●A | |
| PVARP after PVC | 400 ms | Pace | Unipolar |
| AV Search + | Off | Sense | Bipolar |
| Blanking | | ■V | |
| A-Blank after V-Pace | 125 ms | Pace | Unipolar |
| A-Blank after V-Sense | 45 ms | Sense | Bipolar |
| V-Blank after A-Pace | 65 ms | Rate Adaptive Pacing | |
| Magnet Response | Pace Async | Minute Ventilation | Passive |
| Noise Response | DOO | Accelerometer | Passive |
| Rate Enhancements | | | |
| Rate Smoothing | | | |
| Up | Off % | | |
| Down | Off % | | |
| Rate Hysteresis | | | |
| Hysteresis Offset | -5 min ⁻¹ | | |
| Search Hysteresis | Off cycles | | |
| Sudden Brady Response | Off | | |

Setup

Telemetry

Enable use of ZIP™ telemetry On

Trending Setup

Recording Method 30 Second Average
Duration 25 hours
Data Storage Continuous

Sleep Schedule

Sleep Start Time 23:00 hh:mm
Sleep Duration 07 hours

Leads Status Setup

| | Intrinsic Amplitude | Pace Impedance | Pace Impedance Alert Range | Pace Threshold |
|----|---------------------|----------------|----------------------------|----------------|
| ●A | On | On | 200 - 2000 Ω | Off |
| ■V | On | On | 200 - 2000 Ω | Off |