

ZOOM ® View™	Report Created 11 Jan 2016

Combined Follow-up Report

Date of Birth N/R N/R N/R

Device VITALIO MRI J275/101893

Last Office Interrogation

26 Jan 2015 Implant Date N/R

My Alerts

There are no alerts to display.

Events Since Last Reset (02 Jun 2014)

11 Jan 2016 10:57 Tachy at 175 min⁻¹
11 Jan 2016 10:56 Tachy at 166 min⁻¹
See last page for full list of Events Since Last Reset.

Battery OK

Approximate time to explant: 9.5 years
Magnet Rate 100 min⁻¹

One Year Remaining

Explant

Leads Data	Implant N/R N/R N/R	Previous Session	Most Recent
Intrinsic Amplitude	N/R mV	22.9 mV	>25.0 mV @ 158 min ⁻¹
Pace Impedance	N/R Ω	552 Ω	588 Ω
Pace Threshold	N/R V @ N/R ms	Test Failed	1.1 V @ 0.4 ms

Settings

Tachy Settings

Tachy EGM Storage On
Detection Rate 160 min⁻¹

Brady Settings

Mode VVIR Pacing Output

Lower Rate Limit 55 min⁻¹
Maximum Sensor Rate 115 min⁻¹ Sensitivity

Refractory 250 - 250 ms Fixed 2.5 mV

Leads Configuration (Pace/Sense)

Unipolar/Bipolar


2.6 V @ 0.4 ms

Rate Adaptive Pacing

Minute Ventilation On Accelerometer Passive


2869 Software Version: 2.02 J275 Firmware Version: D_v1.02.00(1.01)

Clinician Signature:




Tachy Counters	Since Last Reset 02 Jun 2014	Device Totals	
Episode Counters			
Total Episodes	40.8K	40.8K	
Tachy Episodes	801	801	
Nonsustained Episodes	40.0K	40.0K	
Brady Counters	Reset Before Last N/R day(s) N/R to 02 Jun 2014	Since Last Reset 588 day(s) 02 Jun 2014 to Today	
Counters		·	
% Paced	0	51	
Intrinsic Promotion			
Rate Hysteresis			
% Successful	0	0	

Histograms Reset Before Last (N/R)



Histograms Since Last Reset (02 Jun 2014)



All Events Since Last Reset (02 Jun 2014)

LVCIII3 OIIICC Last Itcs	Ct (02 0dil 2017)
11 Jan 2016 11:03	NonSust at 164 min ⁻¹
11 Jan 2016 11:02	NonSust at 160 min ⁻¹
11 Jan 2016 11:01	NonSust at 173 min ⁻¹
11 Jan 2016 11:00	NonSust at 141 min ⁻¹
11 Jan 2016 10:59	NonSust at 165 min ⁻¹
11 Jan 2016 10:59	NonSust at 178 min ⁻¹
11 Jan 2016 10:59	NonSust at 140 min ⁻¹
11 Jan 2016 10:58	NonSust at 142 min ⁻¹
11 Jan 2016 10:57	Tachy at 175 min⁻¹
11 Jan 2016 10:57	NonSust at 161 min ⁻¹
11 Jan 2016 10:56	Tachy at 166 min⁻¹
11 Jan 2016 10:54	NonSust at 176 min ⁻¹
11 Jan 2016 10:53	NonSust at 159 min ⁻¹
11 Jan 2016 10:53	NonSust at 164 min ⁻¹
11 Jan 2016 10:52	NonSust at 164 min ⁻¹
10 Jan 2016 20:46	Tachy at 131 min⁻¹
10 Jan 2016 20:46	Tachy at 153 min⁻¹
10 Jan 2016 20:38	Tachy at 190 min⁻¹
10 Jan 2016 20:37	Tachy at 170 min⁻¹
10 Jan 2016 20:35	Tachy at 162 min⁻¹
10 Jan 2016 20:34	Tachy at 175 min⁻¹
10 Jan 2016 20:33	Tachy at 163 min⁻¹
10 Jan 2016 20:31	Tachy at 168 min⁻¹
10 Jan 2016 20:31	Tachy at 147 min⁻¹

2869 Software Version: 2.02 J275 Firmware Version: D_v1.02.00(1.01) © 2014

Boston Scientific Corporation
or its affiliates. All rights reserved.

Page 4 of 5

ſ	Clinician Signature:
l	
l	

Events Since Last Reset	(Continued) (02 Jun
10 Jan 2016 20:28	Tachy at 178 min ⁻¹
10 Jan 2016 20:27	Tachy at 185 min⁻¹
10 Jan 2016 20:26	Tachy at 166 min⁻¹
10 Jan 2016 20:23	Tachy at 163 min ⁻¹
10 Jan 2016 20:20	Tachy at 166 min ⁻¹
10 Jan 2016 20:19	Tachy at 152 min ⁻¹
10 Jan 2016 17:07	Tachy at 169 min ⁻¹
10 Jan 2016 17:04	Tachy at 189 min ⁻¹
09 Jan 2016 10:37	Tachy at 168 min⁻¹
09 Jan 2016 09:23	Tachy at 183 min⁻¹
09 Jan 2016 09:22	Tachy at 184 min⁻¹
09 Jan 2016 09:19	Tachy at 168 min ⁻¹
09 Jan 2016 09:18	Tachy at 169 min⁻¹
09 Jan 2016 09:13	Tachy at 174 min⁻¹
09 Jan 2016 09:12	Tachy at 175 min⁻¹
09 Jan 2016 08:50	Tachy at 136 min⁻¹
09 Jan 2016 08:47	Tachy at 175 min⁻¹
09 Jan 2016 08:45	Tachy at 173 min⁻¹
09 Jan 2016 03:43	Tachy at 174 min ⁻¹
27 Dec 2015 17:53	Tachy at 155 min ⁻¹
27 Dec 2015 11:27	Tachy at 186 min ⁻¹
27 Dec 2015 11:25	Tachy at 172 min ⁻¹
27 Dec 2015 11:23	Tachy at 162 min ⁻¹
27 Dec 2015 11:22	Tachy at 186 min ⁻¹
27 Dec 2015 11:18	Tachy at 175 min ⁻¹
27 Dec 2015 11:14	Tachy at 151 min ⁻¹
27 Dec 2015 10:49	Tachy at 183 min ⁻¹
26 Dec 2015 15:46	Tachy at 174 min ⁻¹
26 Dec 2015 12:19	Tachy at 183 min⁻¹
26 Dec 2015 12:09	Tachy at 171 min ⁻¹
26 Dec 2015 12:03	Tachy at 171 min ⁻¹
26 Dec 2015 12:02	Tachy at 169 min ⁻¹
26 Dec 2015 11:59	Tachy at 189 min ⁻¹
26 Dec 2015 11:56	Tachy at 178 min ⁻¹
26 Dec 2015 11:55	Tachy at 175 min ⁻¹
26 Dec 2015 11:52	Tachy at 186 min ⁻¹
26 Dec 2015 11:52	Tachy at 146 min ⁻¹
26 Dec 2015 11:48	Tachy at 168 min ⁻¹
26 Dec 2015 11:43	Tachy at 186 min ⁻¹
26 Jan 2015 13:34	RV Auto

Clinician Signature: