

LATITUDE™ Programming System		
Sensor Trending and Replay Report		
Krajca 540702, Miroslav		

Date of Birth N/R N/R N/R Device

CHARISMA CRT-D G324/ 500899 Tachy Mode Monitor + Therapy


Report Created 23 Feb 2022

Last Office Interrogation 22 Feb 2022 Implant Date

22 Feb 2022

Sensor Trending and Replay Settings	Present	Replay
Lower Rate Limit	75	75 min ⁻¹
Max Sensor Rate	130	130 min⁻¹
Accelerometer	On	On
Response Factor	8	8
Activity Threshold	Medium	Medium
Reaction Time	30	30 s
Recovery Time	2	2 min

Sensor Trending - 30 Second Average



Actual Rate 77 min⁻¹ (Paced) Sensor Replay — — 75 min⁻¹

3868 Software Version: 1.10 G324 Firmware Version: H_v1.00.00(1.01)

© 2014-2019 Boston Scientific Corporation or its affiliates. All rights reserved. Page 1 of 1

Clinician Signature: