

|  | ZOOM ® View™  Marker Legend Report |                          | Report Created 07 Jul 2021 |
|--|------------------------------------|--------------------------|----------------------------|
|  |                                    |                          |                            |
|  |                                    |                          | Last Office Interrogation  |
|  | Date of Birth                      | N/R N/R N/R              | N/R                        |
|  | Device                             | ACCOLADE MRI L310/836133 | Implant Date               |

N/R

AS Atrial Sense - After Refractory
AS-Hy Atrial Sense - In Hysteresis Offset
[AS] Atrial Sense - During Blanking
AP Atrial Pace - Lower Rate

AP↓ Atrial Pace - Rate Smoothing Down
AP↑ Atrial Pace - Rate Smoothing Up
AP-Hy Atrial Pace - At Hysteresis Rate
AP-Sr Atrial Pace - Sensor Rate

AP-Ns Atrial Pace - Noise (asynchronous pacing)

AF Atrial Tachy Zone Sense

VS Ventricular Sense - After Refractory
[VS] Ventricular Sense - Noise First Trigger
VS-Hy Ventricular Sense - At Hysteresis Rate

VP Ventricular Pace - Lower Rate

VP↓ Ventricular Pace - Rate Smoothing Down
VP↑ Ventricular Pace - Rate Smoothing Up
VP-BP Ventricular Pace - Backup Pace
VP-Hy Ventricular Pace - At Hysteresis Rate
VP-Sr Ventricular Pace - Sensor Rate

VP-Ns Ventricular Pace - Noise (asynchronous pacing)
VP-VR Ventricular Pace - Ventricular Rate Regulation

VT VT Zone Sense
AN Atrial Rate Noise
VN Ventricular Rate Noise
Epsd Tachy Start Episode
EpsdEnd Tachy End Episode
Fusion Fusion Beat

LOC Loss of Capture
-- Unclassified Event
#.#V Amplitude Threshold Test
#.#mV Intrinsic Amplitude Test
#.#ms Pulse Width Threshold Test

#### Interval in ms

#### $\Omega$  Lead Impedance Test

2869 Software Version: 2.06 L310 Firmware Version: F\_v1.00.00(1.10) © 2014

Boston Scientific Corporation
or its affiliates. All rights reserved.

Page 1 of 1

Clinician Signature: