

	<b>ZOOM ® View™</b>		<b>Report Created 04 Jun 2014</b>
	<b>Device Settings Report</b>		Last Office Interrogation <b>N/R</b>
	Date of Birth	N/R N/R N/R	Implant Date <b>N/R</b>
	Device	VITALIO MRI J275/101893	

### Programming

Last Programmed 04 Jun 2014

### Tachy

Tachy EGM Storage On  
Detection Rate 160 min<sup>-1</sup> (375 ms)

### Atrial Tachy

#### Therapy

#### Ventricular Regulation

Vent Rate Regulation Off

## Brady

### Settings

Mode VVIR  
Lower Rate Limit 55 min<sup>-1</sup>  
Maximum Sensor Rate 115 min<sup>-1</sup>  
Refractory 250 - 250 ms  
Magnet Response Pace Async  
Noise Response VOO

### Rate Enhancements

Rate Smoothing  
Up Off %  
Down Off %

### Output

Auto 1.1 V @ 0.4 ms

### Sensitivity

Fixed 2.5 mV

### Leads

Pace Unipolar  
Sense Bipolar

### Rate Adaptive Pacing

Minute Ventilation On  
Response Factor 8  
Fitness Level Active  
Ventilatory Threshold 120 min<sup>-1</sup>  
Ventilatory Thresh. Response 70 %  
Accelerometer Passive

## Setup

### Telemetry

Enable use of ZIP™ telemetry On

### Trending Setup

Recording Method 30 Second Average  
Duration 25 hours  
Data Storage Continuous

### Sleep Schedule

Sleep Start Time 23:00 hh:mm  
Sleep Duration 07 hours

### Leads Status Setup

Intrinsic Amplitude	Pace Impedance	Pace Impedance Alert Range	Pace Threshold
On	On	200 - 2000 Ω	On