

LATITUDE™ Programming System Report Created 17 Sep 2021 Marker Legend Report Last Office Interrogation Date of Birth N/R N/R N/R 17 Sep 2021

Implant Date N/R

PROPONENT MRI EL L231/938901

AS Atrial Sense - After Refractory and AFR window

Device

AS-Hy Atrial Sense - In Hysteresis Offset
AS-Fl Atrial Sense - In AFR window
(AS) Atrial Sense - During TARP
[AS] Atrial Sense - During Blanking
AP Atrial Pace - Lower Rate

AP↓ Atrial Pace - Rate Smoothing Down
AP↑ Atrial Pace - Rate Smoothing Up
AP-FB Atrial Pace - Fallback (in ATR)
AP-Hy Atrial Pace - At Hysteresis Rate
AP-Sr Atrial Pace - Sensor Rate
AP→ Atrial Pace - Inserted after AFR

AP-Ns Atrial Pace - Noise (asynchronous pacing)
AP-PP Atrial Pace - Atrial Pacing Preference

AP-PAC Atrial Pace - ProACt

AP-SBR Atrial Pace - Sudden Brady Response
AP-VR Atrial Pace - Ventricular Rate Regulation

PAC PAC

(PAC) PAC in TARP

VS Ventricular Sense - After Refractory
[VS] Ventricular Sense - Noise First Trigger
VS-Hy Ventricular Sense - At Hysteresis Rate
VP Ventricular Pace - Lower Rate or Atrial

VP Ventricular Pace - Lower Rate or Atrial Tracked
VP↓ Ventricular Pace - Rate Smoothing Down
VP↑ Ventricular Pace - Rate Smoothing Up
VP-BP Ventricular Pace - Backup Pace
VP-FB Ventricular Pace - Fallback (in ATR)
VP-Hy Ventricular Pace - At Hysteresis Rate
VP-Sr Ventricular Pace - Sensor Rate

VP-MT Ventricular Pace - Atrial Tracked at MTR
VP-Ns Ventricular Pace - Noise (asynchronous pacing)
VP-VR Ventricular Pace - Ventricular Rate Regulation
VP-SBR Ventricular Pace - Sudden Brady Response

VP-PAC Ventricular Pace - ProACt

VP-PP Ventricular Pace - Atrial Pacing Preference

PVC PVC after Refractory
VT VT Zone Sense
AN Atrial Rate Noise
VN Ventricular Rate Noise

ATR↓ Atrial Tachycardia Sense - Count Down
ATR↑ Atrial Tachycardia Sense - Count Up

ATR-Dur ATR Duration Started
ATR-FB ATR Fallback Started
ATR-End ATR Fallback Ended
FB ATR in Progress
PVP→ PVARP after PVC
PMT-B PMT Termination

V-Epsd Ventricular Tachy Start Episode
V-EpsdEnd Ventricular Tachy End Episode

BTR Brady Tachy Response

RYTHMIQ™ Start/End Episode

Fusion Fusion Beat

© 2014-2017

Boston Scientific Corporation
or its affiliates. All rights reserved.

Page 1 of 2

Clinician Signature:

LOC Loss of Capture
-- Unclassified Event

#.#V Amplitude Threshold Test #.#mV Intrinsic Amplitude Test #.##ms Pulse Width Threshold Test #### Interval (A-A or V-V) in ms #### Ω Lead Impedance Test

3869 Software Version: 1.05 L231 Firmware Version: F_v1.00.00(1.10) Clinician Signature: