

|--|

Combined Follow-up Report

Date of Birth N/R N/R N/R

Device VITALIO MRI J276/102818

Last Office Interrogation

N/R Implant Date N/R

My Alerts

There are no alerts to display.

Events Since Last Reset (29 May 2014)

01 Jun 2014 02:11 ATR at 81 min⁻¹, Avg V Rate in ATR: 84 min⁻¹ 01 Jun 2014 02:06 ATR at 80 min⁻¹, Avg V Rate in ATR: 79 min⁻¹

See last page for full list of Events Since Last Reset.

Battery OK

Approximate time to explant: > 8 years
Magnet Rate 100 min⁻¹

One Year Remaining

Explant

		Explant —————		
Leads Data	Implant N/R N/R N/R	Previous Session	Most Recent	
Atrial				
Intrinsic Amplitude	N/R mV	N/R mV	4.3 mV	
Pace Impedance	N/R Ω	N/R Ω	445 Ω	
Pace Threshold	N/R V @ N/R ms	N/R	Auto 0.8 V @ 0.4 ms	
Ventricular				
Intrinsic Amplitude	N/R mV	N/R mV	>25.0 mV @ 65 min ⁻¹	
Pace Impedance	N/R Ω	N/R Ω	506 Ω	
Pace Threshold	N/R V @ N/R ms	N/R	Auto 0.5 V @ 0.4 ms	
Settings				
Ventricular Tachy Settings	,	Atrial Tachy Settings		

Settings		
Ventricular Tachy Settings		Atrial Tachy Settings
Ventricular Tachy EGM Storage	On	ATR Mode Switch 170 min ⁻¹ DDI
Detection Rate	160 min ⁻¹	
Brady Settings		

Drady Octilligs			
Mode	DDD	Pacing Output	
RYTHMIQ™	Off	Atrial	2.5 V @ 0.4 ms
Lower Rate Limit	55 min⁻¹	Ventricular	2.5 V @ 0.4 ms
Maximum Tracking Rate	130 min ⁻¹	Sensitivity	
Paced AV Delay	80 - 180 ms	Atrial	Fixed 0.75 mV
Sensed AV Delay	65 - 150 ms	Ventricular	Fixed 2.5 mV
A-Refractory (PVARP)	240 - 280 ms	Leads Configuration (Pace/Sense)	
V-Refractory (VRP)	230 - 250 ms	Atrial	Unipolar/Bipolar

Ventricular

Rate Adaptive Pacing
Minute Ventilation Passive
Accelerometer Passive

Unipolar/Bipolar

2869 Software Version: 1.09


J276 Firmware Version: D_v1.02.00(1.01)


© 2010

Boston Scientific Corporation
or its affiliates. All rights reserved.

Page 1 of 6


Clinician Signature:







Ventricular Tachy Counters	Since Last Reset	Device Totals	
Tominous ruony Countries	29 May 2014	201100 101010	
Ventricular Episode Counters	,		
Total Episodes	0	0	
VT Episodes (V>A)	0	0	
SVT Episodes (V<=A)	0	0	
Nonsustained Episodes	0	0	
Brady Counters	Reset Before Last	Since Last Reset	
	N/R	29 May 2014	
Counters			
% A Paced	0	8	
% V Paced	0	99	
Intrinsic Promotion			
AV Search +			
% Successful	0	0	
Rate Hysteresis			
% Successful	0	0	
Atrial Burden			
Atrial Burden %	N/R	<1	
Episodes by Duration			
< 1 minute	0	6	
1 min - < 1 hr	0	0	
1 hr - < 24 hr	0	0	
24 hr - < 48 hr	0	0	
> 48 hr	0	0	
Total PACs	0	19.6K	
Ventricular Counters			
Total PVCs	0	3.1K	
Three or More PVCs	0	0	


Histograms Reset Before Last (N/R)





Histograms Since Last Reset (29 May 2014)





2869 Software Version: 1.09 J276 Firmware Version: D_v1.02.00(1.01) © 2010

Boston Scientific Corporation
or its affiliates. All rights reserved.

Page 5 of 6

Clinician Signature:

All Events Since Last Reset (29 May 2014)

01 Jun 2014 02:11	ATR at 81 min ⁻¹ , Avg V Rate in ATR: 84 min ⁻¹
01 Jun 2014 02:06	ATR at 80 min ⁻¹ , Avg V Rate in ATR: 79 min ⁻¹
31 May 2014 20:35	ATR at 339 min ⁻¹ , Avg V Rate in ATR: 81 min ⁻¹
30 May 2014 18:13	ATR at 95 min ⁻¹ , Avg V Rate in ATR: 71 min ⁻¹
29 May 2014 09:48	ATR at 77 min ⁻¹ , Avg V Rate in ATR: 72 min ⁻¹

2869 Software Version: 1.09 J276 Firmware Version: D_v1.02.00(1.01)

Clinician Signature: