



ZOOM® View™

## Device Settings Report

Report Created 02 Oct 2015

Date of Birth N/R N/R N/R  
Device VITALIO MRI J277/433364

Last Office Interrogation  
N/R  
Implant Date  
N/R

## Programming

Last Programmed 02 Oct 2015

## Ventricular Tachy

Ventricular Tachy EGM Storage On  
Detection Rate 160 min<sup>-1</sup> (375 ms)

## Atrial Tachy

## Therapy

## ATR Mode Switch Details

ATR Mode Switch On  
Trigger Rate 170 min<sup>-1</sup>  
Duration 8 cycles  
Entry Count 8 cycles  
Exit Count 8 cycles  
Fallback  
Mode DDI  
Time 00:30 mm:ss  
ATR Fallback LRL 70 min<sup>-1</sup>  
Ventricular Rate Regulation On  
Maximum Pacing Rate 130 min<sup>-1</sup>

## Atrial Tachy Response

Atrial Flutter Response On  
Trigger Rate 170 min<sup>-1</sup>  
PMT Termination On  
**Atrial Overdrive**  
APP/ProACt Off

## Brady

### Settings

Mode	DDDR	Output	
RYTHMIQ™	Off	●A	Trend 3.5 V @ 0.4 ms
Lower Rate Limit	55 min <sup>-1</sup>	■V	Trend 3.5 V @ 0.4 ms
Maximum Tracking Rate	130 min <sup>-1</sup>	Sensitivity	
Maximum Sensor Rate	130 min <sup>-1</sup>	●A	Fixed 0.75 mV
Paced AV Delay	80 - 180 ms	■V	Fixed 2.5 mV
Sensed AV Delay	65 - 150 ms	Leads	
A-Refractory (PVARP)	240 - 280 ms	●A	
V-Refractory (VRP)	230 - 250 ms	Pace	Bipolar
PVARP after PVC	400 ms	Sense	Bipolar
AV Search +	Off	Safety Switch	On
Blanking		■V	
A-Blank after V-Pace	125 ms	Pace	Bipolar
A-Blank after V-Sense	45 ms	Sense	Bipolar
V-Blank after A-Pace	65 ms	Safety Switch	On
Magnet Response	Pace Async	Rate Adaptive Pacing	
Noise Response	DOO	Minute Ventilation	On
Rate Enhancements		Response Factor	8
Rate Smoothing		Fitness Level	Active
Up	Off %	Ventilatory Threshold	120 min <sup>-1</sup>
Down	Off %	Ventilatory Thresh. Response	70 %
Sudden Brady Response	Off	Accelerometer	Passive

## Setup

### Telemetry

Enable use of ZIP™ telemetry On

### Trending Setup

Recording Method 30 Second Average  
Duration 25 hours  
Data Storage Continuous

### Sleep Schedule

Sleep Start Time 23:00 hh:mm  
Sleep Duration 07 hours

### Leads Status Setup

		Daily Intrinsic Amplitude	Daily Impedance	Impedance Limits Low High	Pace Threshold
●A	Pace/Sense	On	On	200 - 2000 Ω	On
■V	Pace/Sense	On	On	200 - 2000 Ω	On