

ZOOM® View™		Report Created 24 Apr 2018
System Fault Warning		Last Office Interrogation 07 Jun 2017
Date of Birth	N/R N/R N/R	Implant Date
Device	INCEPTA ICD F160/101210	N/R
Tachy Mode	Off	

Boston Scientific

Warning Details

Voltage is too low for projected remaining capacity.

Contact Technical Services with code 1003.

1-800-CARDIAC (North America) or contact your local Boston Scientific representative.

2868 Software Version: 4.04

F160 Firmware

Version: B_v1.02.00(4.01)

© 2016

Boston Scientific Corporation

or its affiliates. All rights reserved.

Page 1 of 1

Clinician Signature: