

ZOOM ® View™	Report Created 27 Jun 2010

**Combined Follow-up Report** 

KRÍZ, Karel

30 Aug 1954

Date of Birth AUTOĞEN X4 CRT-D G179/136339 Device Monitor + Therapy Tachy Mode

Last Office Interrogation

03 Jun 2016 Implant Date 12 Oct 2015

My Alerts

**2**4 Jun 2016 08:54 Atrial Tachycardia Response (ATR) detected for >48 hours.

**27 Jun 2016 11:25** Check LV Lead

**Events Since Last Reset (02 Jun 2016)** 

VF at 268 min<sup>-1</sup>, No Therapy 04 Jun 2016 19:32 VF at 262 min<sup>-1</sup>, No Therapy 04 Jun 2016 19:21 See last page for full list of Events Since Last Reset.

**Battery** OK

> Approximate time to explant: 8 years Charge Time 9.5 s Last Capacitor Re-form

**One Year Remaining** 20 Jan 2016 10:08

**Explant** 

**Implant Previous Session Most Recent Leads Data** 12 Oct 2015 Atrial Intrinsic Amplitude N/R mV 5.1 mV 2.8 mV Pace Impedance  $N/R \Omega$ 786 Ω 855 Ω Pace Threshold N/R V @ N/R ms 1.0 V @ 0.4 ms No data collected Right Ventricular Intrinsic Amplitude N/R mV 10.3 mV 10.1 mV @ 50 min<sup>-1</sup> Pace Impedance  $N/R \Omega$  $355 \Omega$ 424 Ω Pace Threshold N/R V @ N/R ms 1.1 V @ 0.4 ms 1.1 V @ 0.4 ms Left Ventricular N/R mV 4.2 mV Intrinsic Amplitude 7.5 mV Pace Impedance  $N/R \Omega$  $603 \Omega$  $705 \Omega$ Pace Threshold N/R V @ N/R ms 1.1 V @ 0.4 ms 1.0 V @ 0.4 ms **Shock Vector**  $N/R \Omega$ 55 Ω Shock Impedance  $53 \Omega$ Settings

Octungo	
Ventricular Tachy Settings	Atrial Tachy Settings
	001 441 4410 ATDM 1 0 11 4-0

230 min<sup>-1</sup> ATP ATR Mode Switch 170 min<sup>-1</sup> DDI VF 36J, 41J, 41Jx6 VT 185 min<sup>-1</sup> Burst Ramp 31J, 36J, 41Jx4

Brady/CRT Settings			
Mode	DDDR	Pacing Output	
Lower Rate Limit	85 min <sup>−1</sup>	Atrial	2.5 V @ 0.4 ms
Maximum Tracking Rate	130 min <sup>-1</sup>	Right Ventricular	3.0 V @ 0.4 ms
Maximum Sensor Rate	130 min <sup>−1</sup>	Left Ventricular	3.5 V @ 0.4 ms
Paced AV Delay	220 - 250 ms	Sensitivity	
Sensed AV Delay	165 - 190 ms	Atrial	AGC 0.25 mV
A-Refractory (PVARP)	240 - 280 ms	Right Ventricular	AGC 0.6 mV
RV-Refractory (RVRP)	230 - 250 ms	Left Ventricular	AGC 1.0 mV
LV-Refractory (LVRP)	250 ms	Leads Configuration (Pace/Sense)	
Ventricular Pacing Chamber	BiV	Atrial	Bipolar
LV Offset	-60 ms	Right Ventricular	Bipolar
		LV (Pace)	LVRing3>>Can
		LV (Sense)	LVTip1>>LVRing2

2868 Software Version: 3.07

G179 Firmware Version: E\_v1.02.00(2.01)


© 2014 **Boston Scientific Corporation** or its affiliates. All rights reserved. Page 1 of 9


Settings (Continued)
Brady/CRT Settings

Rate Adaptive Pacing Minute Ventilation Accelerometer

On Passive

2868 Software Version: 3.07 G179 Firmware Version: E\_v1.02.00(2.01)






2868 Software Version: 3.07 G179 Firmware Version: E\_v1.02.00(2.01) © 2014


Boston Scientific Corporation
or its affiliates. All rights reserved.


Page 4 of 9




Ventricular Tachy Counters	Since Last Reset 02 Jun 2016	Device Totals	
Ventricular Episode Counters			
Total Episodes	2	6	
Treated			
VF Therapy	0	0	
VT Therapy	0	0	
VT-1 Therapy	0	0	
Commanded Therapy	0	0	
Non-Treated			
No Therapy Programmed	0	0	
Nonsustained Episodes	0	3	
Other Untreated Episodes	2	3	
Ventricular Therapy Counters			
ATP Delivered	0	0	
ATP % Successful	0	0	
Shocks Delivered	0	0	
First Shock % Successful	0	0	
Shocks Diverted	1	1	
Brady/CRT Counters	Reset Before Last	Since Last Reset	
	05 May 2016	02 Jun 2016	
Counters			
% A Paced	28	2	
% RV Paced	94	78	
% LV Paced	97	95	
Intrinsic Promotion			
Rate Hysteresis			
% Successful	0	0	
Atrial Burden			
Atrial Burden %	2	90	
Episodes by Duration			
< 1 minute	1	8	
1 min - < 1 hr	1	3	
1 hr - < 24 hr	3	2	
24 hr - < 48 hr	0	0	
> 48 hr	0	2	
Total PACs	47.9K	10.5K	
Ventricular Counters			
Total PVCs	9.4K	5.5K	
Three or More PVCs	2	0	
Clinical Counters	Since Last Reset 02 Jun 2016	Device Totals	
MRI Counters	0	0	


## **Histograms Reset Before Last (05 May 2016)**





### Histograms Since Last Reset (02 Jun 2016)





2868 Software Version: 3.07 G179 Firmware Version: E\_v1.02.00(2.01) © 2014

Boston Scientific Corporation
or its affiliates. All rights reserved.


Page 7 of 9

Clinician Signature:		

## **Heart Rate Variability**


# Last Measure for 24-hour period ending 28 May 2016 14:09

% of Time Used	70	Mode	DDD
Footprint	18 %	Sensed AV Delay	165 - 190 ms
SDANN	43 ms	Ventricular Pacing Chamber	BiV
		LV Offset	-60 ms



### Reference-measured for 24-hour period ending 13 Oct 2015 14:05

% of Time Used	91	Mode	DDD
Footprint	35 %	Sensed AV Delay	160 - 180 ms
SDANN	118 ms	Ventricular Pacing Chamber	BiV
		LV Offset	0 ms



2868 Software Version: 3.07

G179 Firmware Version: E\_v1.02.00(2.01)

© 2014 Boston Scientific Corporation or its affiliates. All rights reserved. Page 8 of 9

# All Events Since Last Reset (02 Jun 2016)

27 Jun 2016 11:23	ATR at 194 min <sup>-1</sup> , Avg V Rate in ATR: 0 min <sup>-1</sup>
27 Jun 2016 11:23	ATR at 192 min <sup>-1</sup> , Avg V Rate in ATR: 95 min <sup>-1</sup>
27 Jun 2016 11:22	ATR at 194 min <sup>-1</sup> , Avg V Rate in ATR: 95 min <sup>-1</sup>
27 Jun 2016 10:07	ATR at 197 min <sup>-1</sup> , Avg V Rate in ATR: 76 min <sup>-1</sup>
27 Jun 2016 09:13	ATR at 194 min <sup>-1</sup> , Avg V Rate in ATR: 79 min <sup>-1</sup>
27 Jun 2016 09:13	ATR at 192 min <sup>-1</sup> , Avg V Rate in ATR: 95 min <sup>-1</sup>
24 Jun 2016 08:54	ATR at 194 min <sup>-1</sup> , Avg V Rate in ATR: 72 min <sup>-1</sup>
24 Jun 2016 07:38	ATR at 188 min <sup>-1</sup> , Avg V Rate in ATR: 85 min <sup>-1</sup>
24 Jun 2016 07:35	ATR at 186 min <sup>-1</sup> , Avg V Rate in ATR: 93 min <sup>-1</sup>
23 Jun 2016 08:27	ATR at 204 min <sup>-1</sup> , Avg V Rate in ATR: 72 min <sup>-1</sup>
23 Jun 2016 08:26	ATR at 202 min <sup>-1</sup> , Avg V Rate in ATR: 100 min <sup>-1</sup>
23 Jun 2016 08:22	ATR at 190 min <sup>-1</sup> , Avg V Rate in ATR: 87 min <sup>-1</sup>
23 Jun 2016 08:22	ATR at 189 min <sup>-1</sup> , Avg V Rate in ATR: 93 min <sup>-1</sup>
23 Jun 2016 08:21	ATR at 189 min <sup>-1</sup> , Avg V Rate in ATR: 94 min <sup>-1</sup>
05 Jun 2016 18:56	PMT at 130 min <sup>-1</sup>
05 Jun 2016 18:51	PMT at 130 min <sup>-1</sup>
05 Jun 2016 18:48	PMT at 130 min <sup>-1</sup>
05 Jun 2016 18:45	PMT at 130 min <sup>-1</sup>
05 Jun 2016 18:45	PMT at 130 min <sup>-1</sup>
04 Jun 2016 19:32	VF at 268 min <sup>-1</sup> , No Therapy
04 Jun 2016 19:21	VF at 262 min <sup>-1</sup> , No Therapy

2868 Software Version: 3.07 G179 Firmware Version: E\_v1.02.00(2.01) © 2014

Boston Scientific Corporation
or its affiliates. All rights reserved.

Page 9 of 9