



ZOOM ® View™

Device Settings Report

Report Created 02 Oct 2014

Date of Birth N/R N/R N/R
 Device AUTOGEN ICD D174/101582
 Tachy Mode Monitor + Therapy

Last Office Interrogation
N/R
 Implant Date
N/R

Programming

Last Programmed 02 Oct 2014
 Ventricular Tachy Mode 02 Oct 2014 Changed to Monitor + Therapy
 02 Oct 2014 Changed to Monitor Only

Ventricular Tachy

VF 220 min⁻¹ (273 ms)

Detection/Redetection

Initial Duration 1.0 s
 Redetection Dur 1.0 s
 Post-shock Dur 1.0 s

Therapy

QUICK CONVERT™ ATP 250 min⁻¹
 Shock 1 41 J
 Shock 2 41 J
 Additional 41 J Shocks 6

VT 160 min⁻¹ (375 ms)

Detection/Redetection

Initial Duration 10.0 s
 Redetection Dur 1.0 s
 Post-shock Dur 1.0 s

ATP1

Burst 3
 Number of Bursts
 Pulses per Burst
 Initial 10
 Increment 0
 Coupling Interval 81 %
 Decrement 0 ms
 Burst Cycle Length 81 %
 Ramp Decrement 0 ms
 Scan Decrement 0 ms
 Minimum Interval 220 ms

Enhancements

Rhythm ID
 VT Detection On
 Initial Detection On
 Sustained Rate Duration Off mm:ss

Post-Shock Detection

Rhythm ID Setup

Passive Method On
 Active Method On
 Temporary LRL 40 min⁻¹

Common Parameters

Stability (For Post-Shock only) 30 ms
 RhythmMatch™ Threshold 94 %

ATP2

Ramp 3
 Number of Bursts
 Pulses per Burst
 Initial 10
 Increment 0
 Coupling Interval 81 %
 Decrement 0 ms
 Burst Cycle Length 81 %
 Ramp Decrement 10 ms
 Scan Decrement 0 ms
 Minimum Interval 220 ms
 ATP Time-out Off mm:ss

Shocks

Shock 1 31 J
 Shock 2 41 J
 Shock 3 -6 41 J

Ventricular Tachy Therapy Setup

ATP

Ventricular ATP Amplitude 5.0 V
 Ventricular ATP Pulse Width 1.0 ms

Magnet and Beeper

Magnet Response Inhibit Therapy
 Beep During Capacitor Charge Off

Shock (All Shocks)

Waveform Biphasic
 Committed Shock Off
 Lead Polarity Initial
 Shock Lead Vector RV Coil to Can

Atrial Tachy

Therapy

VTR Fallback LRL

70 min⁻¹

Ventricular Regulation

Vent Rate Regulation

Off

Brady

Normal Settings

Mode	VVI	Output	3.5 V @ 0.4 ms
Lower Rate Limit	40 min ⁻¹	■V	
V-Refractory (VRP)	250 ms	Sensitivity	
Noise Response	VOO	■V	AGC 0.6 mV
Rate Enhancements		Leads	
Rate Smoothing		■V	
Up	Off %	Pace	Bipolar
Down	Off %	Sense	Bipolar
Rate Hysteresis		Rate Adaptive Pacing	
Hysteresis Offset	-10 min ⁻¹	Minute Ventilation	Passive
Search Hysteresis	Off cycles	Accelerometer	Passive

Brady (Post-Therapy)

Brady Settings

Lower Rate Limit	60 min ⁻¹	Post Therapy	
Output		Post Therapy Period	00:30 mm:ss
■V	5.0 V @ 1.0 ms		

Setup

Beeper

Beep when Explant is Indicated	On
--------------------------------	----

Telemetry

Enable use of ZIP™ telemetry	On
------------------------------	----

Trending Setup

Recording Method	30 Second Average
Duration	25 hours
Data Storage	Continuous

Sleep Schedule

Sleep Start Time	23:00 hh:mm
Sleep Duration	07 hours

Leads Status Setup

	Daily Intrinsic Amplitude	Daily Impedance	Impedance Limits Low	Impedance Limits High	Beep When Out-of-Range	Pace Threshold
■V Pace/Sense	On	On	200 - 2000 Ω		Off	Off
Shock		On	20 - 125 Ω		Off	